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Enzyme Calibrator

Introduction

- (1) The Enzyme Calibrator is a secondary enzyme reference material (a secondary ERM) for Japan Society of Clinical Chemistry transferable reagents (JSCC transferable reagents). By using this product, routine test values of catalytic concentrations in human serum enzymes can be standardized to the values by the JSCC consensus methods.
- (2) The Enzyme Calibrator is a lyophilized product, which contains human enzymes. The base is bovine serum albumin (BSA). The reactivity of enzymes in this product is equivalent to that in human serum.
- (3) The Enzyme Calibrator carries comparable physical and chemical properties (specific gravity, viscosity, and pH) to those of human serum. Therefore, the pipetting accuracy for this product is comparable to that for human serum.
- (4) The assigned values of AST, ALT, ALP (JSCC Transferable Method, except IFCC Method), CK, LD and γ -GT are determined by using Reference Standard-JSCC Enzyme as a reference measurement standard. The assigned value of ChE is determined by using Reference Standard-ChE for pseudo-ChE as a reference measurement standard. The assigned value of AMY is determined by the JSCC·IFCC transferable method and BG5P substrate method. The assigned value of P-AMY is determined by BG5P substrate method.

Reagent

Enzyme Calibrator 4 bottles × for 3 mL
 Contains 0.004% (when reconstituted) sodium azide.
 Store at 8°C or lower.

Preparation

Accurately add 3.0 mL of deionized or distilled water with a volumetric pipet in a vial of Enzyme Calibrator (for 3 mL) and slowly swirl the vial to completely dissolve the lyophilized cake. The reconstituted solution is stable for 3 days at 2 – 8°C.

Precautions

< Precautions on assays >

- (1) Do not use the reagent described above in any procedures other than those described herein. Refer to package inserts in other reagents used with this product.
- (2) Do not use the reagents described above for any purposes other than those described herein.
- (3) The assigned values in the table below are effective only for the lot number stated herein.

< Precautions for protection from hazards >

- (1) If the reagent comes in contact with mouth, eye or skin, wash off immediately with a large amount of water. Consult a physician if necessary.
- (2) Human blood material was used in the manufacture of Enzyme Calibrator. The raw human material used was tested and found negative for HBsAg, anti-HIV1 & 2 antibodies and anti-HCV antibody. Because no test can offer complete assurance that products derived from human blood will not transmit infectious agents, it is recommended that the reagent be handled with the same biohazard precautions used for patient specimens.
- (3) Be careful not to cut yourself with the aluminum cap when removing from the vial.

< Precautions on disposal >

- (1) When discarding the reagent, dispose of it according to local or national regulations. This product contains zinc chloride (6.7 mg/L as zinc when reconstituted).
- (2) This product contains 0.004% (when reconstituted) of sodium azide as a preservative. Sodium azide may react with copper or lead plumbing to form explosive compounds. Even though this product contains minute quantities of sodium azide, drains should be flushed well with a large amount of water when discarding the solution.

Assigned values

Lot No. **KF741**

Enzyme	Origin	Isozyme	Method for determining the assigned values (Manufacturer's standing measurement procedure)	Assigned values (Catalytic concentration, U/L)
AST	Recombinant *2	S type	JSCC Transferable Method	154 *1
ALT	Recombinant *2	S type	JSCC Transferable Method	148 *1
ALP	Recombinant *2	Liver type	JSCC Transferable Method	448 *1
CK	Recombinant *3	MM type	JSCC Transferable Method	436 *1
ChE	Human plasma	pseudo	JSCC Transferable Method	111
LD	Recombinant *4	I type	JSCC Transferable Method	405 *1
γ -GT	Recombinant *2	II type	JSCC Transferable Method	146 *1
AMY	Human saliva, Recombinant *5	S type, P type	JSCC · IFCC Transferable Method	300 *1
P-AMY	Recombinant *5	P type	Proposed JSCC Transferable Method *6	123
ALP	Recombinant *2	Liver type	IFCC Method	159 *8
AMY	Human saliva, Recombinant *5	S type, P type	BG5P Substrate Method *7	301
P-AMY	Recombinant *5	P type	Immunoinhibition, BG5P Substrate Method *7	137

- *1 AST, ALT, ALP, CK, LD, γ -GT, AMY are Certified under a quality system that conforms with ISO 17034 : 2016.
- *2 Human liver gene
- *3 Human skeletal muscle gene
- *4 Human 1 type gene
- *5 Human pancreas gene
- *6 JSCC : Japanese Journal of Clinical Chemistry, 37 Supp. 1, 25-27 (2008). (in Japanese)
- *7 Calculated by the molar extinction coefficient of 4-nitrophenol.
- *8 To be used as assigned value for L-Type ALP IFCC only.

Ordering information

Code No.	Product	Package	Storage
416-57191	Enzyme Calibrator	4 × for 3 mL	8°C or lower

FUJIFILM Wako Pure Chemical Corporation
 1-2, Doshomachi 3-Chome Chuo-ku, Osaka 540-8605, Japan
 Tel : + 81-6-6203-3749
 Fax : + 81-6-6203-1917
<http://www.wako-chem.co.jp>

Wako

23.03.01K15