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## Immuno-Control 1

## Intended use

- (1) Immuno-Control 1 is the liquid type control serum based on human pooled sera and designed for an accuracy control for CRP, RF or ASO.
- (2) The assigned values are measured by using Wako's reagents for automated analyzers, stated as Assay values.

## Contents

Immuno-Control 1 4 bottles × 2 mL  
Contains human serum.

## Preparation and use

Use this control as supplied.

## Warnings and precautions

### Precautions for protection from hazards

- (1) If the control comes in contact with the mouth, eyes or skin, wash off immediately with a large amount of water. Consult a physician if necessary.
- (2) This product is prepared from human sera which were found to be negative for HBsAg, anti-HIV1 & 2 antibodies and anti-HCV antibody. Because no test can offer complete assurance that product derived from human sera will not transmit infectious agents, it is recommended that this product should be handled with the same precautions used for patient specimens.

### Precautions for use

- (1) Use this control according to the package insert herein.  
Read package inserts of the reagents to be used with this control.
- (2) Store the control under the specified conditions. Do not use control past the expiration date stated on the container label.
- (3) Do not use the control described above for any purpose other than described herein.
- (4) The Assay values of each item is applicable for this lot only.

### Precautions for disposal

- (1) When discarding the control, dispose of them according to local or national regulations.
- (2) The control contains 0.09% sodium azide as a preservative. Sodium azide may react with copper or lead plumbing to form explosive compounds. Even though the control contains minute quantities of sodium azide, drains should be flushed well with a large amount of water, when discarding the control.

### Storage condition and shelf life

- (1) Storage condition : Store at 2-10°C .  
(2) Shelf life : The expiration date stated on the label.

## Ordering information

| Code No.  | Product          | Package  |
|-----------|------------------|----------|
| 417-53101 | Immuno-Control 1 | 4 × 2 mL |

### Assay values

Lot No. KJ027

| Item | Reagent      | Method            | Primary standard           | Unit  | Assay values of Immuno-Control 1 |   |             |
|------|--------------|-------------------|----------------------------|-------|----------------------------------|---|-------------|
| CRP  | LT CRP-HS II | Latex immunoassay | IRMM(ERM-DA470)            | mg/dL | <b>0.82</b>                      | ± | <b>0.08</b> |
| RF   | RF-LT        | Latex immunoassay | NIBSC                      | IU/mL | <b>19</b>                        | ± | <b>2</b>    |
| ASO  | ASO-LT       | Latex immunoassay | Japanese domestic standard | IU/mL | <b>83</b>                        | ± | <b>12</b>   |