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Multi Calibrator Lipids

Intended use

The Wako Multi Calibrator Lipids is to be used with Wako's assay reagents for the determination of total cholesterol (T-Cho), free cholesterol (F-Cho), HDL-cholesterol (HDL-C), LDL-cholesterol (LDL-C), triglyceride (TG) and phospholipids (PL) in specimen.

Contents

Multi Calibrator Lipids 4 bottles × for 3 mL

Contains human serum.

The assigned values are lot dependent.

Determination of assigned values

As a primary standard and a determination method of assigned values of this product are shown as follows.

Preparation and use

Accurately add 3 mL of deionized or distilled water and swirl the vial to dissolve completely the contents of each of the calibrator.

The reconstituted calibrator is stable for 1 week at 2-10°C.

When the calibrator is used with the reagents of HDL-C or LDL-C, the reconstituted calibrator shall be used within 2 days after reconstitution and in storage at 2-10°C.

Warnings and precautions

Precautions for protection from hazards

- (1) If the calibrators come in contact with the mouth, eyes or skin, wash off immediately with a large amount of water. Consult a physician if necessary.
- (2) This product is prepared from human sera which was found to be negative for HBsAg, anti-HIV1 & HIV2 antibodies and anti-HCV antibody. Because no test can offer complete assurance that products derived from human sera will not transmit infectious agents, it is recommended that this product should be handled with the same precautions used for patient specimens.
- (3) Be careful not to cut yourself with the aluminium cap when removing it from the vial.

Precautions for use

- (1) Use this calibrator according to the package insert herein. Read package inserts of the reagents to be used with this calibrator.
- (2) Store the calibrators under the specified conditions. Do not use calibrator past the expiration date stated on the container label.
- (3) Do not use the calibrators described above for any purpose other than described herein.
- (4) The assigned value of each item is lot dependent.

Precautions for disposal

- (1) When discarding the calibrators, dispose of them according to local or national regulations.

Storage condition and shelf life

- (1) Storage condition : Store at 2-10°C.
- (2) Shelf life : The expiration date stated on the label.

Ordering information

| Code No. | Product | Package |
|-----------|--|--------------|
| 464-01601 | Multi Calibrator Lipids (The assigned values are lot dependent) | 4 × for 3 mL |

Assigned values

Lot No. EH962

| Item | Abbreviation | Method | Traceability | Assigned value (mg/dL) |
|-------------------|--------------|---|-------------------------|------------------------|
| Total cholesterol | T-Cho | Enzymatic method | ReCCS (JCCRM 211) | 192 |
| Free cholesterol | F-Cho | Enzymatic method | NIST SRM911 cholesterol | 45.0 |
| HDL-cholesterol | HDL-C | Selective elimination method | ReCCS (JCCRM 224) | 52.5 *2 |
| LDL-cholesterol | LDL-C | Selective elimination method | CDC Reference Method | 116 *2 |
| Triglyceride | TG | Enzymatic method (glycerol blanking method) | ReCCS (JCCRM 224)*1 | 95 |
| Phospholipids | PL | Enzymatic method | Choline Chloride | 214 |

* 1 The value of triglyceride is assigned by the JSCC consensus method.

* 2 When the calibrator is used with the reagents of HDL-C or LDL-C, the reconstituted calibrator shall be used within 2 days after reconstitution and in storage at 2-10°C.