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## Lipids Control Set

### Intended Use

Lipids Control Set is lyophilized product, which is human serum basis, and designed to be used as a quality control material for total cholesterol, HDL-cholesterol, LDL-cholesterol and triglyceride. Use the Lipids Control Set with each reagent of FUJIFILM Wako Pure Chemical Corporation.

### Reagents and Ingredients

Lipids Control 1 2 × for 2 mL

Lipids Control 2 2 × for 2 mL

(Human serum)

[Assigned values are lot dependent.]

### Preparation and use

Accurately add 2.0 mL of ion exchange water or distilled water to the vial and dissolve completely by mixing well before use. The reconstituted solution is stable for 5 days in storage at 2-10°C .

### Precautions for protection from hazards

- (1) If the controls come in contact with mouth, eyes or skin, wash off immediately with a large amount of water. Consult a physician if necessary.
- (2) The controls are prepared from human sera which were found to be negative for HBsAg, anti-HIV1/HIV2 antibodies and anti-HCV antibody.  
Because no test can offer complete assurance that products derived from human blood will not transmit infectious agents, it is recommended that the controls should be handled with the same precautions used for patient specimens.
- (3) Be careful not to cut yourself with the aluminum cap when removing from the vial.

### Precautions for assays

- (1) Use this control according to the package insert herein. Read package insert of the reagents to be used with this product.
- (2) Store the controls under the specified conditions. Do not use controls past the expiration date stated on each control container label.
- (3) Do not use the controls for any purpose other than described herein.
- (4) The Assigned values of each item are applicable for this lot only.

### Precautions for disposal

- (1) When discarding the controls, dispose of them according to local or national regulations.

### Storage condition and shelf life

- (1) Storage condition: 2-10°C
- (2) Shelf life: The expiration date stated on the label.

### Ordering information

Code No.	Product	Package
467-55801	Lipids Control Set	2 × 2 conc. × for 2 mL

### Assigned values

### Lot No. KF817

Item	Method	Reagent	Unit	Lipids Control 1		Lipids Control 2	
				Mean	Standard deviation (SD)	Mean	Standard deviation (SD)
Total cholesterol	Enzymatic Method	L-Type CHO M	mg/dL	157	5.4	213	6.4
HDL-cholesterol	Selective elimination Method	L-Type HDL-C M (3)	mg/dL	41.5	1.3	56.0	1.9
LDL-cholesterol	Selective elimination Method	L-Type LDL-C M	mg/dL	105	2.8	144	4.7
Triglyceride	Enzymatic Method (Glycerol blanking Method)	L-Type Triglyceride M	mg/dL	90	3.4	122	4.7

Concentration of each component of this product is shown as the actually measured value obtained by using the reagents described in the above table. The assay procedures are based on the assay parameters set by us.